

A large orange circle is partially obscured by a horizontal orange bar that spans across the top of the page. The bar has a slight wavy effect and ends with three vertical orange bars of varying heights on the right side.

IVFtech Z-tables

**IVF**tech  
for the next generation

## Z-table

The IVFtech Z tables are designed for use in the operating theater when doing oocyte retrieval. It has a clean smooth heated surface ideal for any warming blocks and media flasks. The draw can be used for storage or the draw can also be heated for heating the tools.

The IVFtech uniformed heated Z tables comes in a wide range of variations.

## Options

Sizes are available in: 600x700mm and 400x600mm

Options for the draws:

- With or without heated draw in the 600x700mm

Options table top:

- Stainless steel or Himacs® tabel top

## HiMacs or Stainless steel:

The Hi-Macs® table tops provides an excellent alternative to Stainless Steel. The Hi-Macs® surface is stain resistant and easy clean. Hi-Macs® is scratch resistant and doesn't scratch as easy as stainless steel, furthermore the Hi-Macs® surface can be refurbished and repolished. Easy to see the writing on the dish when working on the Hi-Macs® surface.

## Advantages

- Easy to move around
- Keep utensils warm in the draw

## Products to be used with the Z-table

- Warming blocks





## **IVFtech**

Toppevadvej 34-38  
DK - 3660 Stenløse

Phone + 45 39 40 25 65  
Fax + 45 39 40 25 64  
[www.ivftech.dk](http://www.ivftech.dk)  
E-mail: [info@ivftech.dk](mailto:info@ivftech.dk)