

A large orange circle is partially visible on the left side of the page. A horizontal orange bar extends from the circle across the top of the page, ending in three vertical orange bars of decreasing size on the right.

Heated tabletop

**IVF**tech  
for the next generation

## Heated tabletop

The Heated table tops are for clean room or laboratory.

Electrical heated area

High precision temperature control makes sure to reach a high success rate for the embryo. the heating table plate is controlled by high precision PID regulator, the heating areas is spited in small sizes to reach the high precision, each part have a separate high precision temperature sensor to monitor and regulate the specific part. Splitting the heating area to small sizes make the total temperature variance much smaller and impended, placing a cold block on the plate will not affect the rest on the heating area, and it makes sure that it have a fast recovery time.

Options

Sizes are available in: 90, 120, 150 and 180cm

Heights are in: 800-850mm

Options table top:

- Standard heat 450 x 650mm
- Full heated 450 x 830mm
- Extra heat 450 x 200mm
- Extra heat 450 x 400
- Witness preparation flush
- Anti-vib
- IVFtech heated glass Ø45mm
- Table top integrated incubator

Advantages

- Easy to use
- Easy to clean
- Big variety of options

Products to be used with IVFtech sterile cabinet

- Warming blocks
- Pull-out shelf
- Mini incubator





## **IVFtech**

Toppevadvej 34-38  
DK - 3660 Stenløse

Phone + 45 39 40 25 65  
Fax + 45 39 40 25 64  
[www.ivftech.dk](http://www.ivftech.dk)  
E-mail: [info@ivftech.dk](mailto:info@ivftech.dk)