

The logo features a large orange circle on the left. A horizontal orange bar extends from the circle across the top of the page. The text "IVFtech Econica" is written in white on the left side of this bar. To the right of the bar, there are three vertical orange bars of varying heights, creating a decorative graphic element.

IVFtech Econica

IVFtech
for the next generation

Econica 0,6 and 1,2

- 1.) New Strap line “Econica the economically affordable cabinet”
- 2.) Simplistic and Cost Effective without compromising the Air Flow quality and distribution as well as consistent heated surfaces associated with IVFtech
- 3.) Compatible with all Major Microscopes
- 4.) Bright and Quiet <51dB
- 5.) Sloping Front – Ergonomic Comfort
- 6.) Adjustable Height Control (Econica 60)
- 7.) Negative Pressure Plenum for greater safety
- 8.) Simple On/Off Fan and Light
- 9.) Scratch resistant glass front?

Front opening:	Fixed
Downflow:	0,25-0,40m/s
Energy consumption:	400W
HEPA filter EN 1822, H-14:	99.995 %MPPS
Prefilter, quality EU 3:	ASHRAE 83
min/max height adjustable:	600-1000mm
min/max height fixed:	800-850mm
Temp. settings:	Room temp. to 50°C
Air quality:	Class A
Voltage-frequency:	120-60 / 230-50
Fuses:	10A
External dimensions:	D500xW650xH975mm
Internal dimensions:	D495xW640xH460
Weight :	120kg
Shipping volume net:	0,6m ³



IVFtech Econica is a standard cabinet.

Options:

1,2m is available with one standard heated area and preparation for microscope or as a cold cabinet.

IVFtech Econica 1,2m is available with a fixed stand.

Front opening:	Fixed
Downflow:	0,25-0,40m/s
Energy consumption:	450W
HEPA filter EN 1822, H-14:	99.995 %MPPS
Prefilter, quality EU 3:	ASHRAE 83
min/max height adjustable:	N/A
min/max height fixed:	800-850mm
Temp. settings:	Room temp. to 50°C
Air quality:	Class A
Voltage-frequency:	120-60 / 230-50
Fuses:	10A
External dimensions:	D635xW1240xH1050mm
Internal dimensions:	D600xW1230XH520mm
Weight :	150kg
Shipping volume net:	1,5m ³





IVFtech

Klintehøj Vænge 3
DK - 3460 Birkerød

Phone + 45 39 40 25 65
Fax + 45 39 40 25 64
www.ivftech.dk
E-mail: info@ivftech.dk